



Mizra Medical

Helping People in Challenging Situations

Personal Male Urinals

Mizra Medical offers a selection of high quality male urinals with many innovative features which allow the user to tailor a urinal to his personal needs. These lightweight systems sold in complete kits. Personal Male Urinals can replace diapers or male external catheters in offering a discrete and comfortable system for controlling incontinence while allowing the user a high degree of mobility.

1. Comfort Urinal:




- Our Comfort Urinal (**DB-770**) is a personal male urinal system which is suited to all users. A delicate diaphragm of latex on the inside of the sheath gently holds the penis in place making **it especially effective for users with small or receding penis**, which may slip out of other personal urinals. The complete kit includes 3 replacement diaphragms, drainage tubes and reusable leg bag.

Mizra Medical Ltd.

21, Hasadna St. Industrial Zone Holon 58814 Israel Telephone: 972-3-5506645 Fax: 972-3-5561105
<mailto:info@mizramedical.com> web: www.mizramedical.com



2. Drip Dry Urinal:



- Our Drip Dry Urinal (**DB-760**) is offered with a harness which is mainly latex. This sturdy personal urinal system stands up to frequent washes and dries very quickly. The complete kit includes drainage tubes and reusable leg bag.

3. Uniflex Urinal:



- Our UniFlex (**DB-750**) Personal Male Urinal has an elastic harness with a variety of Velcro closing tabs. These tabs enable the user to close the harness in different positions to offer both a universal fit and extra ease in removal of the harness. The sheath is attached with snaps and can be detached from the harness for easy cleaning. The complete kit includes drainage tubes and reusable leg bag.